

ATCC medium: 2126 QPX medium

Basal Solution (see below).....900.0 ml
Fetal calf serum (HI).....100.0 ml

Mix and dispense aseptically.

Basal Solution:

MEM (alpha modification, Sigma 0644).....5.1 g
CaCl₂ . 2H₂O1.82 g
KCl.....0.68 g
MgCl₂ . 6H₂O4.36 g
MgSO₄ . 7H₂O3.16 g
NaCl.....24.26 g
HEPES.....5.8 g
Glucose.....0.5 g
Distilled water.....900.0 ml

Adjust pH to 7.2. Filter-sterilize.