

**ATCC medium: 2164 *Haloarcula quadrata* medium**

NaCl.....206.0 g  
MgSO<sub>4</sub> . 7H<sub>2</sub>O .....36.0 g  
KCl.....0.37 g  
CaCl<sub>2</sub> . 2H<sub>2</sub>O .....0.5 g  
MnCl<sub>2</sub>.....0.013 mg  
Yeast extract.....5.0 g  
Agar (if desired).....20.0 g  
Distilled water.....1.0 L

Autoclave at 121C for 15 minutes.

If agar is used, it may be necessary to autoclave the base medium and the salts separately.