

ATCC Medium 2174: Brain Heart Infusion Medium w/ 4µg/ml Vancomycin

Broth Medium

BHI Broth (BD 237500)	37.0 g
DI Water.....	950.0 mL
Vancomycin Solution (See below).....	50 mL

Agar Medium

BHI Agar (BD 211065)	52.0 g
DI Water.....	950.0 mL
Vancomycin Solution (See below).....	50 mL

Autoclave media at 121°C and cool to 50°C

Aseptically add vancomycin solution:


Vancomycin Solution

Vancomycin.....	0.004 g
DI Water.....	50.0 mL

Filter sterilize.

*If the reagent antibiotic powder indicates an assay value of potency, use that value to correctly calculate the amount of powder needed.

Final Product Description



AGAR

BROTH

- Medium should be translucent and tan in color.
- There should not be a precipitate.
- Recommended storage temperature: 2-8°C.