

**ATCC medium: 2218 Potato dextrose agar, quarter-strength**

Infusion from potatoes.....125.0 ml  
Dextrose.....5.0 g  
Agar.....15.0 g  
Distilled water.....875.0 ml

Boil 300.0 g of finely diced potatoes in 500.0 ml of water for 1 hour. Filter through cheese cloth to remove large particulates. Add 125.0 ml of the resulting potato infusion to 875.0 ml of distilled water. Add agar and bring to a boil. Add dextrose. Stir to dissolve dextrose. Dispense as needed. Autoclave at 121C for 15 minutes.