

**ATCC Medium: 2258 R2A Medium**

Yeast Extract.....	0.5 g
Proteose Peptone No. 3 (BD 211693).....	0.5 g
Casamino Acids.....	0.5 g
Glucose.....	0.5 g
Soluble Starch.....	0.5 g
Sodium Pyruvate.....	0.3 g
K <sub>2</sub> HPO <sub>4</sub> .....	0.3 g
MgSO <sub>4</sub> x 7H <sub>2</sub> O.....	0.05 g
Agar.....	15.0 g
DI Water.....	1000 ml

Adjust final pH to 7.2 with potassium phosphate before adding the agar.  
Autoclave at 121°C.