

ATCC medium: 259 Brain heart infusion (BHI) with agar, NaCl, inactivated horse serum, yeast extract, and penicillin

Brain Heart Infusion (BD 237500).....	37.0 g
Agar.....	12.0 g
NaCl.....	20.0 g
Horse serum (inactivated).....	100.0 ml
Yeast Extract (Oxoid 121).....	5.0 g
Penicillin.....	10000 U/ml
Distilled water.....	900.0 ml

Prepare Brain Heart Infusion (BD 237500) according to manufacturer's instructions in 900 ml of water; add agar, NaCl, and yeast extract. Autoclave at 121C for 15 minutes. Aseptically add filter-sterilized horse serum and penicillin to sterile, cooled basal medium.