

ATCC Medium: 2797 YPD w/ 30 mg/L leucine and 30 mg/L uracil

YPD Broth(BD 242820).....	50 g
Leucine.....	30 mg
Uracil.....	30 mg
DI Water.....	1000 ml

Dispense into appropriate vessel and autoclave at 121°C.

For solid medium: Replace YPD Broth with 65 g of YPD Agar (BD 242720) and boil to dissolve before autoclaving.

*YPD is commercially available from BD. To make medium from scratch, follow formulation below:

Yeast Extract.....	10.0 g
Peptone.....	20.0 g
Dextrose.....	20.0 g
Agar (if needed).....	20.0 g
DI Water.....	1000 ml

Adjust for final pH of 5.6. Autoclave at 121°C.