

## **ATCC Medium: 28 Sabouraud Dextrose Agar, Emmons Modification**

### **Agar Medium**

Sabouraud Agar Modified (BD 274720).....50.0 g  
DI Water.....1000 ml

Autoclave at 121° C.

### **Broth Medium**

Sabouraud Broth (BD cat 238230).....30.0g  
DI Water.....1000 ml

Autoclave at 121° C.

\*This product is commercially available from BD.

To make medium from scratch, follow formulation below:

### **Sabouraud Dextrose Agar, Emmons Modification Composition**

Neopeptone.....10.0 g  
Dextrose.....20.0 g  
Agar.....20.0 g  
DI Water.....1000 ml

Final pH 7.0 +/- 0.2.

Omit agar for broth medium.

\*This medium is a modification of regular Sabouraud Dextrose Agar.  
It contains a reduced amount of dextrose and has a neutral pH.