

ATCC medium: 295 S8 medium for *Thiobacillus*

Na₂HPO₄1.2 g
KH₂PO₄1.8 g
MgSO₄ . 7H₂O0.1 g
(NH₄)₂SO₄.....0.1 g
CaCl₂0.03 g
FeCl₃0.02 g
MnSO₄0.02 g
Na₂S₂O₃.....10.0 g
NaHCO₃0.5 g
KNO₃5.0 g
Agar (if needed).....15.0 g
Distilled water.....1.0 L

Autoclave at 121C for 15 minutes.