

ATCC Medium: 42 Desulfovibrio Medium

Peptone.....	5.0 g
Beef Extract.....	3.0 g
Yeast Extract.....	0.2 g
MgSO ₄	1.5 g
Na ₂ SO ₄	1.5 g
Fe(NH ₄) ₂ (SO ₄).....	0.1 g
Glucose.....	5.0 g
Agar.....	15.0 g
Tap water.....	1000 ml

Adjust for pH 7.0 +/- 0.2.

Autoclave medium at 121°C and cool under N₂.

Dispense anaerobically.