

ATCC medium: 430 *Leucothrix* medium

ATCC Medium 429 (see below) with MSG replaced with 0.05% each Tryptone (BD 211705) and Yeast extract.

ATCC Medium 429:

NaCl.....	11.7 g
MgCl ₂ . 6H ₂ O	5.35 g
Na ₂ SO ₄	2.0 g
CaCl ₂ . 2H ₂ O	0.75 g
KCl.....	0.35 g
Tris buffer.....	0.5 g
Na ₂ HPO ₄	0.05 g
Sodium glutamate (MSG).....	10.0 g
Distilled water.....	1.0 L

Adjust pH to 7.6.

Autoclave at 121C for 15 minutes.