

ATCC medium: 479 Hagem's-Modess medium

Agar, Noble (BD 214230).....15.0 g
Glucose.....10.0 g
K₂HPO₄0.15 g
KH₂PO₄0.35 g
MgSO₄ . 7H₂O0.50 g
FeCl₃ or Ferric citrate.....1.0 mg
DL-Asparagine.....1.0 g
Yeast extract.....1.0 g
Peptone.....1.0 g
Thiamine.....40.0 mg
Distilled water.....1.0 L

Autoclave at 121C for 15 minutes.