

**ATCC medium: 522 Peptone Czapek's**

NaNO <sub>3</sub> .....	3.0 g
K <sub>2</sub> HPO <sub>4</sub> .....	1.0 g
MgSO <sub>4</sub> .....	0.5 g
KCl .....	0.5 g
FeSO <sub>4</sub> .....	0.01 g
Peptone .....	5.0 g
Sucrose .....	15.0 g
Agar .....	15.0 g
Distilled water .....	1.0 L

Autoclave at 121C for 15 minutes.