

**ATCC medium: 557 Casamino-peptone-Czapek medium**

Casamino Acids (BD 223050).....1.0 g  
Peptone.....2.0 g  
K<sub>2</sub>HPO<sub>4</sub> .....1.0 g  
KCl.....0.5 g  
MgSO<sub>4</sub> . 7H<sub>2</sub>O.....0.5 g  
FeSO<sub>4</sub> . 7H<sub>2</sub>O .....0.01 g  
Sucrose.....30.0 g  
Agar.....15.0 g  
Distilled deionized water.....1.0 L

Autoclave at 121C for 15 minutes.