

ATCC Medium 595: TNT Medium

Agar/Broth Medium

Tryptone.....	10.0 g
NaCl.....	5.0 g
Thiaminehydrochloride.....	1.0 mg
Agar (if required).....	15.0 g
DI Water.....	1.0 L

Stab Medium

Tryptone.....	10.0 g
NaCl.....	5.0 g
Thiaminehydrochloride.....	1.0 mg
Agar (0.5%).....	5.0 g
DI Water.....	1.0 L

Adjust pH of all media to 7.3.

Autoclave at 121°C and dispense as required.