

ATCC medium: 607 Na/suc medium with carbenicillin

Basal medium:

Brain Heart Infusion (BD 237500)	37.0 g
Yeast extract	5.0 g
NaCl	9.3 g
Sucrose	97.3 g
MgSO ₄	0.25 g
Agar (if needed)	13.3 g
Deionized water to	1.0 L

Adjust pH to 7.4. Autoclave at 121C for 15 minutes.

To 900 ml of this solution, aseptically add:

Horse serum, gamma globulin-free (heat-inactivated)	100.0 ml
Disodium carbenicillin	0.5 mg