

**ATCC medium: 609 *Pseudomonas* medium**

Nitrilotriacetic acid.....1.91 g  
K<sub>2</sub>HPO<sub>4</sub> .....8.71 g  
Na<sub>2</sub>SO<sub>4</sub> .....0.57 g  
MgSO<sub>4</sub> .....0.25 g  
FeSO<sub>4</sub> .....0.5 mg  
Ca (NO<sub>3</sub>)<sub>2</sub> .....0.5 mg  
Agar.....15.0 g  
Distilled water.....1.0 L

Adjust pH to 6.5.  
Autoclave at 121C for 15 minutes.