

ATCC medium: 62 *Escherichia* medium

Trypticase Peptone (BD 211921)	6.0 g
K ₂ HPO ₄	0.2 g
FeSO ₄ . 7H ₂ O	5.0 mg
MgSO ₄ . 7H ₂ O	0.2 g
L-Asparagine	0.15 g
Glycerol	2.0 g
Vitamin B12	0.4 mg
Agar (if needed)	15.0 g
Distilled water	1.0 L

Adjust medium for final pH 7.2. Autoclave at 121C for 15 minutes.