

ATCC medium: 698 Peptone-succinate agar (ATCC medium 234) with seawater instead of distilled water

ATCC Medium 234:

(NH ₄) ₂ SO ₄	1.0 g
MgSO ₄ . 7H ₂ O	1.0 g
MnSO ₄ . H ₂ O	2.0 mg
FeCl ₃ . 6H ₂ O	2.0 mg
Succinic acid.....	1.68 g
Peptone.....	5.0 g
Agar.....	1.5 g
Distilled water.....	1.0 L

Adjust final pH to 7.0.