

ATCC medium: 730 *Microbacterium* medium

FeSO ₄ . 7H ₂ O	10.0	mg
Calcium D-(+)-pantothenate.....	200.0	mg
Thiamine . HCl.....	10.0	mg
(NH ₄) ₂ SO ₄	2.0	g
KH ₂ PO ₄	5.0	g
K ₂ HPO ₄	5.0	g
Potassium aspartate.....	5.0	g
Glucose.....	10.0	g
MgSO ₄ . 7H ₂ O	500.0	mg
Biotin.....	0.1	mg
beta-Mercaptopurine.....	100.0	mg
Distilled water.....	1.0	L

Adjust pH to 7.0. Autoclave 10 minutes at 121C.