

ATCC medium: 768 YEPD-FA medium

Dextrose.....	20.0 g
Peptone.....	20.0 g
Yeast extract.....	10.0 g
K ₂ HPO ₄	5.0 g
KH ₂ PO ₄	5.0 g
Tween 40.....	10.0 g
Myristic acid.....	0.07 g
Palmitic acid.....	0.07 g
Stearic acid.....	0.07 g
Agar.....	20.0 g
Distilled water.....	1.0 L

Autoclave at 121C for 15 minutes.