

ATCC medium: 779 G medium

MgSO₄0.2 g
K₂HPO₄0.6 g
KH₂PO₄0.4 g
(NH₄)₂SO₄.....2.0 g
Yeast extract.....2.0 g
MnSO₄ . H₂O0.05 g
ZnSO₄ . 7H₂O5.0 mg
CuSO₄ . 5H₂O5.0 mg
CaCl₂0.08 g
FeSO₄ . 7H₂O0.5 mg
Glucose.....1.0 g
Distilled water.....1.0 L

Adjust for a final pH to 7.8.
Autoclave at 121C for 15 minutes.