

ATCC medium: 788 Benzene sulfonate medium

Benzene sulfonic acid, sodium salt.....	1.0 g
(NH ₄) ₂ SO ₄	1.0 g
MgSO ₄ . 7H ₂ O	0.2 g
KH ₂ PO ₄	0.3 g
K ₂ HPO ₄	0.7 g
CaCl ₂	10.0 mg
FeSO ₄ . 7H ₂ O	5.0 mg
MnSO ₄ . 5H ₂ O	2.0 mcg
CuSO ₄	50.0 mcg
H ₃ BO ₃	10.0 mcg
ZnSO ₄ . 7H ₂ O	70.0 mcg
MoO ₃ . 2H ₂ O	10.0 mcg
Agar.....	15.0 g
Distilled water.....	1.0 L

Autoclave at 121C for 15 minutes.