

**ATCC medium: 78 *Lactobacillus* medium**

Agar.....	15.0 g
Tryptone.....	20.0 g
Tryptose.....	5.0 g
Yeast extract.....	5.0 g
Tomato juice (filtered, pH 7.0).....	200.0 ml
Liver Extract Concentrate (ICN 100377).....	1.0 g
Tween 80.....	0.05 g
Glucose.....	3.0 g
Lactose.....	2.0 g
Distilled water to.....	1.0 L

Adjust pH to 6.5.

Autoclave at 121C for 15 minutes.