

**ATCC Medium 8: Nutrient Agar/Broth w/ 3% NaCl**

**Agar**

Nutrient Agar (BD 213000).....	23.0 g
NaCl.....	30.0 g
DI Water.....	1000.0 ml

**Broth**

Nutrient Broth (BD 234000).....	8.0 g
NaCl.....	30.0 g
DI Water.....	1000.0 ml

Autoclave at 121°C. Dispense as required.

\*Nutrient Agar/Broth is commercially available from BD, to make from scratch, follow formulation below.

**BD Nutrient Agar/Broth Composition**

Beef Extract.....	3.0 g
Peptone.....	5.0 g
Agar (if required).....	15.0 g

Final pH 6.8 +/- 0.2.