

**ATCC Medium: 8 Nutrient Agar/Broth w/ 3% NaCl**

**For Solid Media**

Nutrient Agar (BD 213000).....23.0 g  
NaCl.....30.0 g  
DI Water.....1000 ml

Autoclave at 121°C. Final pH 6.8 +/- 0.2.

**For Broth**

Nutrient Broth (BD 234000).....8.0 g  
NaCl .....30.00 g  
DI Water ..... 1000 mL

Autoclave at 121C. Final pH 6.8 +/- 0.2.