

ATCC medium: 950 Pontecorvo's *Aspergillus* medium (ATCC medium 687) with 0.05% yeast extract

ATCC Medium 687 (see below) with 0.05% yeast extract

ATCC Medium 687:

NaNO ₃	6.0 g
KH ₂ PO ₄	1.52 g
KCl	0.52 g
MgSO ₄ . 7H ₂ O	0.52 g
ZnSO ₄	0.001 g
FeCl ₃	0.85 mg
Glucose	30.0 g
Trace Metals (see below)	1.0 ml
Biotin	40.0 mcg
Distilled water	1.0 L

Adjust pH to 6.2 with KOH.

Trace Metals:

Na ₂ B ₄ O ₇ . 10H ₂ O	100.0 mg
MnCl ₂ . 4H ₂ O	50.0 mg
Na ₂ MoO ₄ . 2H ₂ O	50.0 mg
CuSO ₄ . 5H ₂ O	250.0 mg
Distilled water	1.0 L