

**ATCC medium: 1660 MH salts**

NaCl.....	120.5 g
MgCl <sub>2</sub> . 6H <sub>2</sub> O .....	22.4 g
MgSO <sub>4</sub> .....	14.4 g
CaCl <sub>2</sub> .....	0.54 g
KCl.....	3.0 g
NaHCO <sub>3</sub> .....	0.09 g
NaBr.....	39.0 mg
Proteose Peptone No. 3 (BD 211693).....	5.0 g
Yeast extract.....	10.0 g
Glucose.....	1.0 g
Agar.....	20.0 g
Distilled water to.....	1.0 L

Adjust medium for final pH 7.5. Autoclave at 121C for 15 minutes.