

ATCC Medium: 1705 Brucella Agar/Broth with 5% Defibrinated Sheep Blood

Agar Medium

Brucella Agar (BD 211086).....43.0 g
Sheep Blood (defibrinated; aseptically added).....50 ml
DI Water.....950 ml

Autoclave at 121°C. Cool sterilized medium to ~47°C. Aseptically add 50 ml of room temperature defibrinated sheep blood. Gently mix and dispense as required.

Broth Medium

Brucella Broth (BD 211088).....28.0 g
Sheep Blood (defibrinated; aseptically added).....50 ml
DI Water.....950 ml

Autoclave at 121°C. Cool sterilized medium to ~47°C. Aseptically add 50 ml of room temperature defibrinated sheep blood. Gently mix and dispense as required.