

ATCC medium: 894 TEP Uric acid

Na ₂ HPO ₄ . 12H ₂ O.....	9.0 g
KH ₂ PO ₄	1.5 g
MgSO ₄ . 7H ₂ O	0.2 g
Ferric ammonium citrate.....	1.2 mg
CaCl ₂	20.0 mg
MnCl ₂ . 4H ₂ O	1.0 mg
Uric acid.....	4.0 g
Trypticase Soy Broth (BD 211768).....	3.0 g
Agar.....	20.0 g
Distilled water.....	1.0 L

Adjust pH to 7.2.

Autoclave at 121C for 15 minutes.