


Technical Data Sheet: THP-1 AP1-LUC2

ATCC [®] Number	TIB-202-AP1-LUC2™
Organism	Homo sapiens
Tissue/Disease Source	Acute Monocytic Leukemia
Product Description	THP-1 cell line (ATCC TIB-202 [™]) is commonly used to study human monocyte and macrophage activities, functions, innate immune mechanisms and signaling pathways. This luciferase reporter cell line was derived from parental line TIB-202 by stably expressing firefly luciferase gene (luc2) under control of an activator protein 1 (AP1) promoter. This cell line was established through lentiviral transduction and single cell cloning. The cells, upon stimulation, express high levels of enzymatically active luciferase protein, which can be detected via <i>in vitro</i> bioluminescence assays. This reporter cell line is useful for monitoring the activity of cell signal transduction pathways in immune cells, specifically detecting the activation of the JNK signaling pathway are triggered by stress, inflammatory cytokines, mitogens and oncogenes.
Application	Enabling sensitive and quantitative assessment of signal transduction makes this reporter cell line ideal for <i>in vitro</i> bioluminescence assays to study immune response in human monocytes, development of new drugs, and safety evaluation of new chemicals and drugs.



Luciferase Expression

Figure 1: Linearity of luminescence and in vitro quantification of luciferase activity of THP-1 **AP1-LUC2.** Cells were seeded in a 96-well plate, after stimulating overnight, bioluminescence signals were detected using Bright-Glo[™] (Promega[®]) and a luminometer (Glomax[™]). Error bars show standard deviation (n=3). (A) Cells were seeded at indicated cell numbers and stimulated with 10 µg/mL of PMA to determined the linear correlation of bioluminescence intensity with cell numbers. (B) Serial dilution of stimulus reagents (PMA or LPS) to demonstrate dose response sensitivity of the AP1 promoter.

Cell Morphology

THP-1 (TIB-202™)

THP-1 AP1-LUC2 (TIB-202-AP1-LUC2™)




Figure 2: Cell morphology of THP-1 parental and THP-1 AP1-LUC2. Cells were maintained in ATCC recommended culture conditions. Cell morphology was observed under microscopy and images were captured by digital camera. Red bar represents 1000 µm.

© 2020 American Type Culture Collection. The ATCC trademark and trade name, and any other trademarks listed in this publication are owned by the American Type Culture Collection unless indicated otherwise.

These products are for laboratory use only. Not for human or diagnostic use. ATCC products may not be resold, modified for resale, used to provide commercial services, or to manufacture commercial products without prior ATCC written approval.