



lambdafoo

87099™

Product Sheet

Description

Clone type: Vector

Storage Conditions

Product format: Freeze-dried

Intended Use

This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use.

BSL 1

ATCC determines the biosafety level of a material based on our risk assessment as guided by the current edition of *Biosafety in Microbiological and Biomedical Laboratories (BMBL)*, U.S. Department of Health and Human Services. It is your responsibility to understand the hazards associated with the material per your organization's policies and procedures as well as any other applicable regulations as enforced by your local or national agencies.

Certificate of Analysis

For batch-specific test results, refer to the applicable certificate of analysis that can

be found at www.atcc.org.

Vector Information

Construct size (kb): 43.0

Intact vector size: 43.000

Vector name: lambdafoo (phage, lambda - replacement)

Type of vector: phage

Construction: lambda2001, lambda gene V

Host range: *Escherichia coli*

Vector information:

other: amber stop codon, ->

other: Pro-Thr box encoding a linker peptide, ->

other: amber stop codon

other: Pro-Thr box encoding a linker peptide

Cloning sites: HindIII; BamHI; SacI; EcoRI

Coding sequence: gene V N-terminal 176 aa, -> gene V N-terminal 176 aa

Initiation codon: ATG

Insert detection: lacZ', ->

MCS: HindIII...EcoRI, ->

Polylinker sites: HindIII; BamHI; SacI; EcoRI

Restriction sites: SfiI; SfiI

Ribosome-binding site: Shine-Dalgarno sequence

Growth Conditions

Medium:

ATCC Medium 1592: SM buffer

Temperature: 37°C

Notes

Restriction digests of the clone give the following sizes (kb): BamHI--33.0, 9.4; EcoRI--33.0, 9.4; BglII--22.0, 8.8, 4.8, 4.6, 3.1; PstI--9.6, 9.0, 4.6,

3.2, 2.9, 2.8, 2.4, 2.2, 1.9, 1.9, 1.5.

- ATCC staff

Vector allowing expression of cloned inserts as a fusion protein on the phage particle surface. Inserts are fused to the C-terminus of a truncated phage tail protein by a peptide linker.

- Proc. Natl. Acad. Sci. USA 91: 8273-8277, 1994

Presence of the lacZalpha coding sequence and ribosome binding site allow blue-white color detection of recombinants, as well as allowing expression of a cloned insert separate from the phage tail protein.

- Proc. Natl. Acad. Sci. USA 91: 8273-8277, 1994

The Pro-Thr box encodes alternating prolyl and threonyl residues, which form a link between the N-terminal phage tail protein and the foreign protein.

- Proc. Natl. Acad. Sci. USA 91: 8273-8277, 1994

The linker resembles the hinge-region of IgA1 and allows separation of the foreign protein from the phage particles by digestion with enzymes such as Cellulomonas fimi protease or collagenase.

- Proc. Natl. Acad. Sci. USA 91: 8273-8277, 1994

Production of large amounts of fusion protein may inhibit phage assembly. A host allowing low efficiency suppression of the amber mutation is recommended.

- Proc. Natl. Acad. Sci. USA 91: 8273-8277, 1994

Material Citation

If use of this material results in a scientific publication, please cite the material in the following manner: lambdafoo (ATCC 87099)

References

References and other information relating to this material are available at www.atcc.org.

Warranty

The product is provided 'AS IS' and the viability of ATCC® products is warranted for 30 days from the date of shipment, provided that the customer has stored and handled the product according to the information included on the product information sheet, website, and Certificate of Analysis. For living cultures, ATCC lists the media formulation and reagents that have been found to be effective for the product. While other unspecified media and reagents may also produce satisfactory results, a change in the ATCC and/or depositor-recommended protocols may affect the recovery, growth, and/or function of the product. If an alternative medium formulation or reagent is used, the ATCC warranty for viability is no longer valid. Except as expressly set forth herein, no other warranties of any kind are provided, express or implied, including, but not limited to, any implied warranties of merchantability, fitness for a particular purpose, manufacture according to cGMP standards, typicality, safety, accuracy, and/or noninfringement.

Disclaimers

This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use. Any proposed commercial use is prohibited without a [license from ATCC](#).

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate or complete and the customer bears the sole responsibility

of confirming the accuracy and completeness of any such information.

This product is sent on the condition that the customer is responsible for and assumes all risk and responsibility in connection with the receipt, handling, storage, disposal, and use of the ATCC product including without limitation taking all appropriate safety and handling precautions to minimize health or environmental risk. As a condition of receiving the material, the customer agrees that any activity undertaken with the ATCC product and any progeny or modifications will be conducted in compliance with all applicable laws, regulations, and guidelines. This product is provided 'AS IS' with no representations or warranties whatsoever except as expressly set forth herein and in no event shall ATCC, its parents, subsidiaries, directors, officers, agents, employees, assigns, successors, and affiliates be liable for indirect, special, incidental, or consequential damages of any kind in connection with or arising out of the customer's use of the product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of such materials.

Please see the material transfer agreement (MTA) for further details regarding the use of this product. The MTA is available at www.atcc.org.

Copyright and Trademark Information

© ATCC 2023. All rights reserved.

ATCC is a registered trademark of the American Type Culture Collection.

Revision

This information on this document was last updated on 2025-09-12

Contact Information

ATCC

10801 University Boulevard

lambdafoo

87099

Manassas, VA 20110-2209

USA

US telephone: 800-638-6597

Worldwide telephone: +1-703-365-2700

Email: tech@atcc.org or contact your local distributor

Product Sheet