



Haemophilus haemoglobinophilus (Neumann) Murray

19416™

Description

- **Strain designation** NCTC 1659 [CIP 53.88 T, XIII]
 - **Deposited As** *Haemophilus canis* (sic) Rivers
 - **Type strain** Yes
-

Storage Conditions

- **Product format** Freeze-dried
-

Intended Use

This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use.

BSL 1

ATCC determines the biosafety level of a material based on our risk assessment as guided by the current edition of *Biosafety in Microbiological and Biomedical Laboratories (BMBL)*, U.S. Department of Health and Human Services. It is your responsibility to understand the hazards associated with the material per your organization's policies and procedures as well as any other applicable regulations as enforced by your local or national agencies.

ATCC highly recommends that appropriate personal protective equipment is always used when handling vials. For cultures that require storage in liquid nitrogen, it is important to note that some vials may leak when submersed in liquid nitrogen and will slowly fill with liquid nitrogen. Upon thawing, the conversion of the liquid nitrogen back to its gas phase may result in the vial exploding or blowing off its cap with dangerous force creating flying debris. Unless necessary, ATCC recommends that these cultures be stored in the vapor phase of liquid nitrogen rather than submersed

Haemophilus haemoglobinophilus (Lehmann and Neumann) Murray

¹⁹⁴¹⁶ in liquid nitrogen.

Certificate of Analysis

For batch-specific test results, refer to the applicable certificate of analysis that can be found at www.atcc.org.

Growth Conditions

- **Medium**
[ATCC Medium 814: GC Agar/Broth Medium](#)
 - **Temperature** 37°C
 - **Atmosphere** 95% Air, 5% CO₂
-

Handling Procedures

- **1. Open vial according to enclosed instructions.**
 - **2. Using 0.5 ml of #18 broth, rehydrate the vial contents.**
 - **3. Aseptically transfer 0.2 ml of this suspension to #814 slant(s), and/or 0.1 ml to plates.**
 - **4. Incubate tubes and plates at 37°C in an atmosphere of 5% CO₂.**
 - **5. Within 3 to 4 days, growth should be evident by colonies on the agar surfaces. Cells are Gram negative, thin regular to filamentous rods. Colonies are small pinpoint circular with a convex elevation, an entire margin, and a smooth surface.**
-

Notes

If a broth culture is desired, use Medium #814 formulation without the agar, or HTM (Haemophilus Test Broth). To make HTM, first a fresh hematin stock solution is prepared by dissolving 50 mg of hematin powder in 100 ml of 0.01N (0.01 mol/L) NaOH with heat and stirring until the powder is dissolved thoroughly. Thirty milliliters of the hematin stock solution is added to

Haemophilus haemoglobinophilus (Lehmann and Neumann) Murray

¹⁹⁴¹⁶ 1 L of Mueller-Hinton Broth with 5 g of yeast extract. After autoclaving and cooling, aseptically add 3 ml of a NAD (nicotinamide adenine dinucleotide) stock solution (50 mg of NAD dissolved in 10 ml of distilled water; filter sterilized).

Material Citation

If use of this material results in a scientific publication, please cite the material in the following manner: *Haemophilus haemoglobinophilus* (Lehmann and Neumann) Murray (ATCC 19416)

References

References and other information relating to this material are available at www.atcc.org.

Warranty

The product is provided 'AS IS' and the viability of ATCC® products is warranted for 30 days from the date of shipment, provided that the customer has stored and handled the product according to the information included on the product information sheet, website, and Certificate of Analysis. For living cultures, ATCC lists the media formulation and reagents that have been found to be effective for the product. While other unspecified media and reagents may also produce satisfactory results, a change in the ATCC and/or depositor-recommended protocols may affect the recovery, growth, and/or function of the product. If an alternative medium formulation or reagent is used, the ATCC warranty for viability is no longer valid. Except as expressly set forth herein, no other warranties of any kind are provided, express or implied, including, but not limited to, any implied warranties of merchantability, fitness for a particular purpose, manufacture according to cGMP standards, typicality, safety, accuracy, and/or noninfringement.

Disclaimers

- This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use. Any proposed commercial use is prohibited without a license from ATCC.

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate or complete and the customer bears the sole

***Haemophilus haemoglobinophilus* (Lehmann and Neumann) Murray**

¹⁹⁴¹⁶ responsibility of confirming the accuracy and completeness of any such information.

This product is sent on the condition that the customer is responsible for and assumes all risk and responsibility in connection with the receipt, handling, storage, disposal, and use of the ATCC product including without limitation taking all appropriate safety and handling precautions to minimize health or environmental risk. As a condition of receiving the material, the customer agrees that any activity undertaken with the ATCC product and any progeny or modifications will be conducted in compliance with all applicable laws, regulations, and guidelines. This product is provided 'AS IS' with no representations or warranties whatsoever except as expressly set forth herein and in no event shall ATCC, its parents, subsidiaries, directors, officers, agents, employees, assigns, successors, and affiliates be liable for indirect, special, incidental, or consequential damages of any kind in connection with or arising out of the customer's use of the product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of such materials.

Please see the material transfer agreement (MTA) for further details regarding the use of this product. The MTA is available at www.atcc.org.

Copyright and Trademark Information

© ATCC 2023. All rights reserved.
ATCC is a registered trademark of the American Type Culture Collection.

Revision

This information on this document was last updated on 2021-05-19

Contact Information

ATCC
10801 University Boulevard
Manassas, VA 20110-2209
USA
US telephone: 800-638-6597
Worldwide telephone: +1-703-365-2700
Email: tech@atcc.org or contact your local distributor