



pJP1 plasmid in *Bacillus sphaericus*

67080™

Description

This clone contains the complete coding sequence (nt 245-1414 of Genbank M15686) of the gene lysostaphin from *Staphylococcus simulans*. It was constructed by excising the insert of pRG5 (ATCC 67076) with an EcoRI/HindIII digest, and blunt-ending the vector and insert with the Klenow fragment of DNA polymerase before ligation.

– U.S. Patent 4,931,390 dated June 5, 1990.

Organism: *Staphylococcus simulans* Kloos and Schleifer

Clone type: Clone

Host: *Bacillus sphaericus*

Patent depository: This material was deposited with the ATCC Patent Depository to fulfill U.S. or international patent requirements. This material may not have been produced or characterized by ATCC. As an International Depository Authority (IDA) for patent deposits, ATCC is required to complete viability testing only at time of initial deposit of patent material. Patent deposits are made available on behalf of the Depositor when the pertinent U.S. or international patent is issued, but material may not be used to infringe the patent claims.

Patent number:

4,931,390

Technical information: ATCC Technical Services does not have technical information on patent deposits that are not produced or characterized by ATCC. Additional information can be found in the corresponding patent available from the patent holder or with the U.S. and/or international patent office.

Shipping information: *B. sphaericus* containing the plasmid

Storage Conditions

Product format: Freeze-dried

Storage conditions: 2°C to 8°C

Intended Use

This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use.

BSL 1

ATCC determines the biosafety level of a material based on our risk assessment as guided by the current edition of *Biosafety in Microbiological and Biomedical Laboratories (BMBL)*, U.S. Department of Health and Human Services. It is your responsibility to understand the hazards associated with the material per your organization's policies and procedures as well as any other applicable regulations as enforced by your local or national agencies.

Certificate of Analysis

For batch-specific test results, refer to the applicable certificate of analysis that can be found at www.atcc.org.

Insert Information

Insert size (kb): 1.5

Insert information:

Gene: lysostaphin

Source: *Staphylococcus simulans*, NRRL B-2628

Insert ends: (5') EcoRI ; (3 ') HindIII

Genbank accession: M15686

Insert end: 3' HindIII; 5' EcoR

Vector Information

pJP1 plasmid in *Bacillus sphaericus* 67080

Construct size (kb): 8.0

Intact vector size: 4.6

Vector name: pBC16

Type of vector: plasmid

Vector end: EcoRI

Markers: eryR

Growth Conditions

Medium:

ATCC Medium 1179: VY medium

Handling Procedures

1. Open vial according to instructions.
 2. Aseptically add 0.3 to 0.4 mL of liquid medium to the freeze-dried pellet and mix well. Transfer 100 μ L to a test tube containing 5 mL VY Medium (*see below) + erythromycin (5 μ g/mL). A loopful of culture can also be streaked on an agar plate of the same. Incubate cultures at 30°C.
 3. Isolate DNA using standard plasmid preparation procedures for *Bacillus* sp.
-

Notes

Restriction digests of the clone gave the following sizes (in kb): EcoRI/HindIII – 8.0 ;

PstI – uncut ; Sph – 8.9.

–ATCC Staff

Material Citation

If use of this material results in a scientific publication, please cite the material in the following manner: pJP1 plasmid in *Bacillus sphaericus* (ATCC 67080)

References

References and other information relating to this material are available at www.atcc.org.

Warranty

The product is provided 'AS IS' and the viability of ATCC® products is warranted for 30 days from the date of shipment, provided that the customer has stored and handled the product according to the information included on the product information sheet, website, and Certificate of Analysis. For living cultures, ATCC lists the media formulation and reagents that have been found to be effective for the product. While other unspecified media and reagents may also produce satisfactory results, a change in the ATCC and/or depositor-recommended protocols may affect the recovery, growth, and/or function of the product. If an alternative medium formulation or reagent is used, the ATCC warranty for viability is no longer valid. Except as expressly set forth herein, no other warranties of any kind are provided, express or implied, including, but not limited to, any implied warranties of merchantability, fitness for a particular purpose, manufacture according to cGMP standards, typicality, safety, accuracy, and/or noninfringement.

Disclaimers

This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use. Any proposed commercial use is prohibited without a license from ATCC.

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate or complete and the customer bears the sole responsibility

of confirming the accuracy and completeness of any such information.

This product is sent on the condition that the customer is responsible for and assumes all risk and responsibility in connection with the receipt, handling, storage, disposal, and use of the ATCC product including without limitation taking all appropriate safety and handling precautions to minimize health or environmental risk. As a condition of receiving the material, the customer agrees that any activity undertaken with the ATCC product and any progeny or modifications will be conducted in compliance with all applicable laws, regulations, and guidelines. This product is provided 'AS IS' with no representations or warranties whatsoever except as expressly set forth herein and in no event shall ATCC, its parents, subsidiaries, directors, officers, agents, employees, assigns, successors, and affiliates be liable for indirect, special, incidental, or consequential damages of any kind in connection with or arising out of the customer's use of the product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of such materials.

Please see the material transfer agreement (MTA) for further details regarding the use of this product. The MTA is available at www.atcc.org.

This material is cited in a US and/or international patent and may not be used to infringe the claims. Depending on the wishes of the Depositor, ATCC may be required to inform the Depositor of the party to which the material was furnished.

Copyright and Trademark Information

© ATCC 2023. All rights reserved.

ATCC is a registered trademark of the American Type Culture Collection.

Revision

This information on this document was last updated on 2022-10-22

Contact Information

pJP1 plasmid in *Bacillus sphaericus*
67080

ATCC

10801 University Boulevard

Manassas, VA 20110-2209

USA

US telephone: 800-638-6597

Worldwide telephone: +1-703-365-2700

Email: tech@atcc.org or contact your local distributor
