



clone 2.72F, #19

95491™

Description

Detects:

Genome: Homo sapiens

Tissue: fibroblast

Gene symbol: ITGA2

Type of nucleic acid: mRNA

Invariant bands (kb): 8.0

Organism: *Homo sapiens*, human

Clone type: Clone

Shipping information: Rehydrate with TE

Storage Conditions

Product format: Frozen

Intended Use

This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use.

BSL 1

ATCC determines the biosafety level of a material based on our risk assessment as guided by the current edition of *Biosafety in Microbiological and Biomedical Laboratories (BMBL)*, U.S. Department of Health and Human Services. It is your responsibility to understand the hazards associated with



clone 2.72F, #19

95491

the material per your organization's policies and procedures as well as any other applicable regulations as enforced by your local or national agencies.

Certificate of Analysis

For batch-specific test results, refer to the applicable certificate of analysis that can be found at www.atcc.org.

Insert Information**Insert size (kb):** 5.367**Type of DNA:** cDNA**Insert source:** lung fibroblast IMR-90 cell line**Insert tissue:** lung fibroblast IMR-90 cell line**Insert information:**

DESCRIPTION OF INSERT COMPONENT:

Genomic copy number: unique
receptor)

Cross references: DNA Seq. Acc.: M28249

Nucleotides 1-5367 of the insert correspond to
nucleotides 2-5368 of M28249.**Genome:** Homo sapiens**Chromosome:** 5

5 q23-q31

Gene name: integrin, alpha 2 (CD49B, alpha 2 subunit of VLA-2**Gene product:** integrin, alpha 2 (CD49B, alpha 2 subunit of VLA-2
receptor)(collagen receptor, platelet gpla/IIa) [ITGA2]**Gene symbol:** ITGA2**Contains complete coding sequence:** Yes**Insert end:** Modification: EcoRI linker**Vector Information**

clone 2.72F, #19

95491

Construct size (kb): 8.329999923706055**Intact vector size:** 2.964**Vector name:** pBluescript KS+**Type of vector:** phagemid**Construction:** pUC19**Host range:** *Escherichia coli***Vector end:** EcoRI closer to XhoI site; EcoRI**Cloning sites:** SacII; XmaI; NotI; XbaI; SpeI; BamHI; SmaI; PstI; EcoRI; EcoRV; HindIII; ClaI; Sall HincII Accl; XhoI; DralI; ApaI; KpnI**Enhancer:** none**Insert detection:** lacZ'**Markers:** ampR**Polylinker sites:** SEE COMMENTS**Promoters:** T7; T3; lac**Replicon:** pMB1; f1**Terminator:** none

Notes

Restriction digests of the clone give the following sizes (kb): EcoRI--3.6, 3.0, 1.8; HindIII--3.2, 3.2, 2.1; XbaI--8.4; XhoI--4.8, 3.5; KpnI--4.8, 3.5.

- ATCC staff

The insert contains the following restriction sites (approximate kb from the 5'

end): EcoRI--1.79, BamHI--1.86; HindIII--2.02, 5.2.

- J. Cell Biol. 109: 397-407, 1989

The coding region is nt 49-3594.

- J. Cell Biol. 109: 397-407, 1989

Material Citation

clone 2.72F, #19

95491

If use of this material results in a scientific publication, please cite the material in the following manner: clone 2.72F, #19 (ATCC 95491)

References

References and other information relating to this material are available at www.atcc.org.

Warranty

The product is provided 'AS IS' and the viability of ATCC® products is warranted for 30 days from the date of shipment, provided that the customer has stored and handled the product according to the information included on the product information sheet, website, and Certificate of Analysis. For living cultures, ATCC lists the media formulation and reagents that have been found to be effective for the product. While other unspecified media and reagents may also produce satisfactory results, a change in the ATCC and/or depositor-recommended protocols may affect the recovery, growth, and/or function of the product. If an alternative medium formulation or reagent is used, the ATCC warranty for viability is no longer valid. Except as expressly set forth herein, no other warranties of any kind are provided, express or implied, including, but not limited to, any implied warranties of merchantability, fitness for a particular purpose, manufacture according to cGMP standards, typicality, safety, accuracy, and/or noninfringement.

Disclaimers

This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use. Any proposed commercial use is

prohibited without a license from ATCC.

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate or complete and the customer bears the sole responsibility of confirming the accuracy and completeness of any such information.

This product is sent on the condition that the customer is responsible for and assumes all risk and responsibility in connection with the receipt, handling, storage, disposal, and use of the ATCC product including without limitation taking all appropriate safety and handling precautions to minimize health or environmental risk. As a condition of receiving the material, the customer agrees that any activity undertaken with the ATCC product and any progeny or modifications will be conducted in compliance with all applicable laws, regulations, and guidelines. This product is provided 'AS IS' with no representations or warranties whatsoever except as expressly set forth herein and in no event shall ATCC, its parents, subsidiaries, directors, officers, agents, employees, assigns, successors, and affiliates be liable for indirect, special, incidental, or consequential damages of any kind in connection with or arising out of the customer's use of the product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of such materials.

Please see the material transfer agreement (MTA) for further details regarding the use of this product. The MTA is available at www.atcc.org.

Copyright and Trademark Information

© ATCC 2021. All rights reserved.

ATCC is a registered trademark of the American Type Culture Collection.

Revision

clone 2.72F, #19

95491

This information on this document was last updated on 2021-05-19

Contact Information

ATCC

10801 University Boulevard

Manassas, VA 20110-2209

USA

US telephone: 800-638-6597

Worldwide telephone: +1-703-365-2700

Fax number: 703-365-2701

Email: tech@atcc.org or contact your local distributor
