



p58-CL49

95706™

Description

Organism: *Homo sapiens*, human

Clone type: Clone

Host: *Escherichia coli* HB101 (ATCC 33694)

Storage Conditions

Product format: Freeze-dried

Intended Use

This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use.



Biosafety Level 1

ATCC determines the biosafety level of a material based on our risk assessment as guided by the current edition of *Biosafety in Microbiological and Biomedical Laboratories (BMBL)*, U.S. Department of Health and Human Services. It is your responsibility to understand the hazards associated with the material per your organization's policies and procedures as well as any other applicable regulations as enforced by your local or national agencies.

Certificate of Analysis

For batch-specific test results, refer to the applicable certificate of analysis that can be found at www.atcc.org.

Insert Information

Insert size (kb): 1.5529999999999999

Type of DNA: cDNA

Insert source: Resting CD3- cells

Insert tissue: Resting CD3- cells

Insert information:

DESCRIPTION OF INSERT COMPONENT:

ORF's seq. position: GDB probe: GDB:698358 ()

Insert 5' end: Modification: Sall linker

Insert 3' end: Modification: NotI linker

Cross references: DNA Seq. Acc.: U24079

Nucleotides 1-1535 of the insert correspond to nucleotides 1-1535 of U24079.

Genome: Homo sapiens

Chromosome: 19

19 q13.4; 19, 374.860000 cR

Gene name: killer cell inhibitory receptor

Gene product: killer cell inhibitory receptor [KIR]

Gene symbol: KIR; NKB1

Contains complete coding sequence: Yes

Vector Information

Construct size (kb): 5.650000095367432

Intact vector size: 4.109

Vector name: pSPORT1

Type of vector: plasmid

Host range: *Escherichia coli*

Vector end: Sall; NotI

Cloning sites: AatII; SphI; mLuI; SunI; SnaBI; HindIII; BamHI; XbaI; XmaIII; NotI; SpeI; SstI; Sall; Accl; SmaI; EcoRI; Kpn2I; RsrII; PinAI; KpnI; Sse8387I; PstI

Insert detection: lacZ', <-, 1-341

Markers: ampR

MCS: AatII...PstI, ->, 191-285

Polylinker sites: AatII; SphI; mLuI; SunI; SnaBI; HindIII; BamHI; XbaI; XmaIII; NotI; SpeI; SstI; Sall; Accl; SmaI; EcoRI; Kpn2I; RsrII; PinAI; KpnI; Sse8387I; PstI

Promoters: SP6; T7

Replicon: pMB1, 1903-1903

Repressor gene: lacI, <-, 464-1546

Growth Conditions

Medium:

ATCC Medium 1227: LB Medium (ATCC medium 1065) with 50 mcg/ml ampicillin

Temperature: 37°C

Notes

Restriction digests of the clone give the following sizes (kb): NotI/Sall--4.4, 1.6; EcoRV--5.7; PstI--4.8, 0.4, 0.35, 0.1; PvuII--3.7, 1.4, 0.6, 0.1; XhoI--5.7. Recommended for verification: Sall, NotI, NotI/Sall, EcoRI, PstI.
- ATCC staff

Encodes a variant recognized by monoclonal antibody GL183 but not EB6.
- Immunity 2: 439-449, 1995

Insert contains the following restriction sites (approximate kb from the 5' end): PstI--0.04, 0.44, 0.51, 0.84; PvuII--1.13; SacI--0.45; StuI--1.06; XhoI--0.63.
- Immunity 2: 439-449, 1995

Material Citation

If use of this material results in a scientific publication, please cite the material in the following manner: p58-CL49 (ATCC 95706)

References

References and other information relating to this material are available at www.atcc.org.

Warranty

The product is provided 'AS IS' and the viability of ATCC® products is warranted for 30 days from the date of shipment, provided that the customer has stored and handled the product according to the information included on the product information sheet, website, and Certificate of Analysis. For living cultures, ATCC lists the media formulation and reagents that have been found to be effective for the product. While other unspecified media and reagents may also produce satisfactory results, a change in the ATCC and/or depositor-recommended protocols may affect the recovery, growth, and/or function of the product. If an alternative medium formulation or reagent is used, the ATCC warranty for viability is no longer valid. Except as expressly set forth herein, no other warranties of any kind are provided, express or implied, including, but not limited to, any implied warranties of merchantability, fitness for a particular purpose, manufacture according to cGMP standards, typicality, safety, accuracy, and/or noninfringement.

Disclaimers

This product is intended for laboratory research use only. It is not intended

for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use. Any proposed commercial use is prohibited without a license from ATCC.

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate or complete and the customer bears the sole responsibility of confirming the accuracy and completeness of any such information.

This product is sent on the condition that the customer is responsible for and assumes all risk and responsibility in connection with the receipt, handling, storage, disposal, and use of the ATCC product including without limitation taking all appropriate safety and handling precautions to minimize health or environmental risk. As a condition of receiving the material, the customer agrees that any activity undertaken with the ATCC product and any progeny or modifications will be conducted in compliance with all applicable laws, regulations, and guidelines. This product is provided 'AS IS' with no representations or warranties whatsoever except as expressly set forth herein and in no event shall ATCC, its parents, subsidiaries, directors, officers, agents, employees, assigns, successors, and affiliates be liable for indirect, special, incidental, or consequential damages of any kind in connection with or arising out of the customer's use of the product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of such materials.

Please see the material transfer agreement (MTA) for further details regarding the use of this product. The MTA is available at www.atcc.org.

Copyright and Trademark Information

© ATCC 2021. All rights reserved.

ATCC is a registered trademark of the American Type Culture Collection.

Revision

This information on this document was last updated on 2021-05-19

Contact Information

ATCC

10801 University Boulevard

Manassas, VA 20110-2209

USA

US telephone: 800-638-6597

Worldwide telephone: +1-703-365-2700

Fax number: 703-365-2701

Email: tech@atcc.org or contact your local distributor
