



# Genomic DNA from *Trichoderma virens* anamorph strain train T-1

9645D-2™

## Description

Genomic DNA isolated from *Trichoderma virens* T-1. This fungal strain is also available as ATCC® Catalog No.: 9645.

**Organism:** *Trichoderma virens* (Miller et al.) von Arx

**Derived from:** *Trichoderma virens* T-1 [CBS 430.54, DSM 1963, IAM 5061, IFO 6355, IMI 45553, NRRL 2314, QM 365] (ATCC 9645)

**Genome sequenced strain:** Yes

**Type strain:** No

**Mass:** 2 µg

**Shipping information:** Stored in 1X TE buffer

---

## Storage Conditions

**Product format:** Freeze-dried

---

## Intended Use

This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use.

---



BSL 1

# Genomic DNA from *Trichoderma virens* anamorph strain T-1

9645D-2

ATCC determines the biosafety level of a material based on our risk assessment as guided by the current edition of Biosafety in Microbiological and Biomedical Laboratories (BMBL), U.S. Department of Health and Human Services. It is your responsibility to understand the hazards associated with the material per your organization's policies and procedures as well as any other applicable regulations as enforced by your local or national agencies.

---

## Certificate of Analysis

For batch-specific test results, refer to the applicable certificate of analysis that can be found at [www.atcc.org](http://www.atcc.org).

---

## Handling Procedures

### Centrifuge tube prior to opening to prevent loss of pelleted material

1. Rehydrate contents of vial with molecular grade H<sub>2</sub>O.  
DNA is dried in 1X Tris buffer.
  2. Place vial at 37°C for 1 hour or at +2°C to 8°C overnight.
  3. For more complete rehydration and to fully recover DNA incubate the sample overnight at 4°C while rocking, then incubate for 1 hour at 65°C.
- 

## Quality Control Specifications

**Electrophoresis - RNA content:** No RNA was detected by electrophoresis

**Purity (A260/A280):** 1.7 to 2.1

**Integrity:** Integrity of DNA was determined by electrophoresis on a 1% agarose gel stained with SYBR Safe™, and was found to be of high molecular weight.

**Functional tests:** Functional activity was confirmed by PCR amplification of approximately 1500 base pairs fragment of rRNA gene cluster including ITS1-

# Genomic DNA from *Trichoderma virens* anamorph strain train T-1

9645D-2

5.8S-ITS2 region.

**Identity:** Identity confirmed by sequencing of ITS1, 5.8S gene and ITS2 regions of ribosomal RNA (~ 500 base pairs).

---

## Notes

Genomic DNA isolated from fungi is appropriate for PCR\* and other molecular biology applications.

\*The polymerase chain reaction (PCR) process is covered by patents owned by Hoffmann-LaRoche Inc. Use of the PCR process requires a license.

---

## Material Citation

If use of this material results in a scientific publication, please cite the material in the following manner: Genomic DNA from *Trichoderma virens* anamorph strain train T-1 (ATCC 9645D-2)

---

## References

References and other information relating to this material are available at [www.atcc.org](http://www.atcc.org).

---

## Warranty

The product is provided 'AS IS' and the viability of ATCC® products is warranted for 30 days from the date of shipment, provided that the customer has stored and handled the product according to the information included on the product information sheet, website, and Certificate of Analysis. For living cultures, ATCC lists the media formulation and reagents that have been found to be effective for the product. While other

9645D-2

unspecified media and reagents may also produce satisfactory results, a change in the ATCC and/or depositor-recommended protocols may affect the recovery, growth, and/or function of the product. If an alternative medium formulation or reagent is used, the ATCC warranty for viability is no longer valid. Except as expressly set forth herein, no other warranties of any kind are provided, express or implied, including, but not limited to, any implied warranties of merchantability, fitness for a particular purpose, manufacture according to cGMP standards, typicality, safety, accuracy, and/or noninfringement.

---

## Disclaimers

This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use. Any proposed commercial use is prohibited without a license from ATCC.

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate or complete and the customer bears the sole responsibility of confirming the accuracy and completeness of any such information.

This product is sent on the condition that the customer is responsible for and assumes all risk and responsibility in connection with the receipt, handling, storage, disposal, and use of the ATCC product including without limitation taking all appropriate safety and handling precautions to minimize health or environmental risk. As a condition of receiving the material, the customer agrees that any activity undertaken with the ATCC product and any progeny or modifications will be conducted in compliance with all applicable laws, regulations, and guidelines. This product is provided 'AS IS' with no representations or warranties whatsoever except as expressly set forth herein and in no event shall ATCC, its parents, subsidiaries, directors, officers,

9645D-2

agents, employees, assigns, successors, and affiliates be liable for indirect, special, incidental, or consequential damages of any kind in connection with or arising out of the customer's use of the product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of such materials.

Please see the material transfer agreement (MTA) for further details regarding the use of this product. The MTA is available at [www.atcc.org](http://www.atcc.org).

---

### **Copyright and Trademark Information**

© ATCC 2021. All rights reserved.

ATCC is a registered trademark of the American Type Culture Collection.

---

### **Revision**

This information on this document was last updated on 2021-05-19

---

### **Contact Information**

ATCC

10801 University Boulevard

Manassas, VA 20110-2209

USA

US telephone: 800-638-6597

Worldwide telephone: +1-703-365-2700

Fax number: 703-365-2701

Email: [tech@atcc.org](mailto:tech@atcc.org) or contact your local distributor

---