



FL 3-4 1.5RI:0.5RI [FL3/41.5RI:0.5RI]

99534™

Description

Organism: *Homo sapiens*, human

Clone type: Clone

Host: *Escherichia coli* HB101 (ATCC 33694)

Storage Conditions

Product format: Freeze-dried

Intended Use

This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use.

BSL 1

ATCC determines the biosafety level of a material based on our risk assessment as guided by the current edition of *Biosafety in Microbiological and Biomedical Laboratories (BMBL)*, U.S. Department of Health and Human Services. It is your responsibility to understand the hazards associated with the material per your organization's policies and procedures as well as any other applicable regulations as enforced by your local or national agencies.

Certificate of Analysis

For batch-specific test results, refer to the applicable certificate of analysis that can be found at www.atcc.org.

Insert Information

Insert size (kb): 1.8

Type of DNA: cDNA

Insert information:

DESCRIPTION OF INSERT COMPONENT:

Insert 5' end: XbaI

Insert 3' end: EcoRI

Cross references: DNA Seq. Acc.: U12707

Nucleotides 1-1783 of the insert correspond to nucleotides 24-1806 of U12707.

Genome: Homo sapiens

Chromosome: X

X p11.3-p11.22

Gene name: Wiskott-Aldrich syndrome (eczema-thrombocytopenia)

Gene product: Wiskott-Aldrich syndrome (eczema-thrombocytopenia)(Wiskott-Aldrich syndrome protein) [WAS]

Gene symbol: WAS; WASP; IMD2

Contains complete coding sequence: Yes

Vector Information

Construct size (kb): 4.760000228881836

Intact vector size: 2.960

Vector name: pBluescript SK-

Type of vector: phagemid

Construction: pUC19

Host range: *Escherichia coli*

Vector end: XbaI; EcoRI

Cloning sites: KpnI; ApaI; DraI; XhoI; Sall; HincII; AclI; ClaI; HindIII; EcoRV; EcoRI; PstI;

FL 3-4 1.5RI:0.5RI [FL3/41.5RI:0.5RI]

99534

SmaI; BamHI; SpeI; XbaI; EagI; NotI; SacII; BstXI; SacI

Insert detection: lacZ', <-**Markers:** ampR**MCS:** KpnI...SacI, ->**Polylinker sites:** SacI; BstXI; SacII; NotI; EagI; XbaI; SpeI; BamHI; SmaI; PstI; EcoRI;

EcoRV; HindIII; (ClaI Bsp106I); (HincII AccI SalI); XhoI; (Eco0109I DralI ApaI); KpnI

Primer site: M13-20, -> M13 reverse, <-**Promoters:** T3; T7; lac**Replicon:** f1, →; pMB1**Growth Conditions****Medium:**

ATCC Medium 1227: LB Medium (ATCC medium 1065) with 50 mcg/ml ampicillin

Temperature: 37°C**Notes**

Restriction digests of the clone give the following sizes (kb): SalI/XbaI--3.0, 1.8; EcoRI--4.4, 0.5; PstI--3.4, 1.2; XhoI--3.6, 0.8.

- ATCC staff

The insert contains the following restriction sites (approximate kb from the 5' end): PstI--0.15, 1.38; BglII--0.81; XhoI--0.97; EcoRI--1.33.

- GenBank/EMBL/DDBJ

Constructed by ligating a 0.5 kb EcoRI fragment corresponding to the 3' end of the insert to an existing construct.

- personal communication

Muscle, liver, pancreas, kidney, testis, lung, heart, brain, eye and prostate do not express WAS.

- Cell 78: 635-644, 1994

Material Citation

If use of this material results in a scientific publication, please cite the material in the following manner: FL 3-4 1.5RI:0.5RI [FL3/41.5RI:0.5RI] (ATCC 99534)

References

References and other information relating to this material are available at www.atcc.org.

Warranty

The product is provided 'AS IS' and the viability of ATCC® products is warranted for 30 days from the date of shipment, provided that the customer has stored and handled the product according to the information included on the product information sheet, website, and Certificate of Analysis. For living cultures, ATCC lists the media formulation and reagents that have been found to be effective for the product. While other unspecified media and reagents may also produce satisfactory results, a change in the ATCC and/or depositor-recommended protocols may affect the recovery, growth, and/or function of the product. If an alternative medium formulation or reagent is used, the ATCC warranty for viability is no longer valid. Except as expressly set forth herein, no other warranties of any kind are provided, express or implied, including, but not limited to, any implied warranties of merchantability, fitness for a particular purpose, manufacture according to cGMP standards, typicality, safety, accuracy, and/or noninfringement.

Disclaimers

99534

This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use. Any proposed commercial use is prohibited without a [license from ATCC](#).

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate or complete and the customer bears the sole responsibility of confirming the accuracy and completeness of any such information.

This product is sent on the condition that the customer is responsible for and assumes all risk and responsibility in connection with the receipt, handling, storage, disposal, and use of the ATCC product including without limitation taking all appropriate safety and handling precautions to minimize health or environmental risk. As a condition of receiving the material, the customer agrees that any activity undertaken with the ATCC product and any progeny or modifications will be conducted in compliance with all applicable laws, regulations, and guidelines. This product is provided 'AS IS' with no representations or warranties whatsoever except as expressly set forth herein and in no event shall ATCC, its parents, subsidiaries, directors, officers, agents, employees, assigns, successors, and affiliates be liable for indirect, special, incidental, or consequential damages of any kind in connection with or arising out of the customer's use of the product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of such materials.

Please see the material transfer agreement (MTA) for further details regarding the use of this product. The MTA is available at www.atcc.org.

Copyright and Trademark Information

© ATCC 2023. All rights reserved.

ATCC is a registered trademark of the American Type Culture Collection.

Revision

FL 3-4 1.5RI:0.5RI [FL3/41.5RI:0.5RI]

99534

Product Sheet

This information on this document was last updated on 2024-10-25

Contact Information

ATCC

10801 University Boulevard

Manassas, VA 20110-2209

USA

US telephone: 800-638-6597

Worldwide telephone: +1-703-365-2700

Email: tech@atcc.org or contact your local distributor
