

PSC-derived CD34+ Cells, BXS0117

CS-7020™

Description

iPSC-derived CD34+ Cells, BXS0117 are induced pluripotent stem cells that have been differentiated into hematopoietic progenitor CD34+ cells. These cells can be further differentiated down common lymphoid progenitor and common myeloid progenitor lineages. iPSC-derived CD34+ Cells can be used in cancer immunology research, drug development, toxicity screening, and blood lineage differentiation studies.

- **Organism** *Homo sapiens*, human
- **Age** 27 years
- **Gender** Female
- **Morphology** rounded
- **Growth properties** Suspension
- **Cells per vial** Approximately 3.0×10^6
- **Volume** 1.0 mL

Storage Conditions

- **Product format** Frozen
- **Storage conditions** Vapor phase of liquid nitrogen

Intended Use

This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use.

BSL 2

ATCC determines the biosafety level of a material based on our risk assessment as guided by the current edition of *Biosafety in Microbiological and Biomedical Laboratories (BMBL)*, U.S. Department of Health and Human Services. It is your responsibility to understand the hazards associated with the material per your organization's policies and procedures as well as any other applicable regulations as enforced by your local or national agencies.

ACS-7020

ATCC highly recommends that appropriate personal protective equipment is always used when handling vials. For cultures that require storage in liquid nitrogen, it is important to note that some vials may leak when submersed in liquid nitrogen and will slowly fill with liquid nitrogen. Upon thawing, the conversion of the liquid nitrogen back to its gas phase may result in the vial exploding or blowing off its cap with dangerous force creating flying debris. Unless necessary, ATCC recommends that these cultures be stored in the vapor phase of liquid nitrogen rather than submersed in liquid nitrogen.

Certificate of Analysis

For batch-specific test results, refer to the applicable certificate of analysis that can be found at www.atcc.org.

Growth Conditions

- **Temperature** 37°C
- **Atmosphere** 95% Air, 5% CO₂

Handling Procedures

- **Unpacking and storage instructions**
 1. Check all containers for leakage or breakage.
 2. Remove the frozen cells from the dry ice packaging and immediately place the cells at a temperature below -130°C, preferably in liquid nitrogen vapor, until ready for use.
- **Complete medium** iPSC-derived CD34+ hematopoietic progenitor cells should be thawed prior to their intended use in application specific media. ATCC recommends thawing them in RPMI-1640 (ATCC 30-2001). ATCC does not recommend maintaining iPSC-derived CD34+ hematopoietic progenitor cells in culture in the absence of application-specific growth factors.
- **Handling Procedure** Refer to the batch specific information for the total number of viable cells recovered from this lot of ATCC ACS-7020 Lot 70020687 is $\sim 2.8 \times 10^6$
 1. Using the total number of viable cells, customers have to decide seeding for their experiments and applications ².
 2. Prepare the desired combinations of culture dishes with required media. Place dishes in a 37°C, 5% CO₂, humidified incubator and allow the media to pre-equilibrate to temperature

ACS-7020

and pH for 30 minutes prior to adding cells.

3. While the culture dishes equilibrate, remove one vial of ATCC ACS-7020 from storage and thaw the cells by gentle agitation in a 37°C water bath. To reduce the possibility of contamination, keep the O-ring and cap out of the water. Thawing should be rapid (approximately 1 to 2 minutes).
 4. Remove the vial from the water bath as soon as the contents are thawed, and decontaminate by dipping in or spraying with 70% ethanol. All operations from this point onward should be carried out under strict aseptic conditions.
 5. Add 4ml of base medium or complete growth media – into a sterile conical tube. Using a sterile pipette, transfer cells from the cryovial to the conical tube. Centrifuge at 200-300xg for 5min, remove supernatant and re suspend the pellet in complete growth medium
 6. Transfer cell suspension to each of the pre-equilibrated culture dishes in the required seeding density, gently rock each dishes to evenly distribute the cells.
 7. Place the seeded culture flasks in the incubator at 37°C with a 5% CO₂ atmosphere. Incubate for at least 24 hours before processing the cells further.
- **Cryopreservation N/A:** As this cell line is intended to be consumable no sub culturing and no cryopreservation is recommended.
-

Material Citation

If use of this material results in a scientific publication, please cite the material in the following manner: iPSC-derived CD34+ Cells, BXS0117 (ATCC ACS-7020)

References

References and other information relating to this material are available at www.atcc.org.

Warranty

The product is provided 'AS IS' and the viability of ATCC® products is warranted for 30 days from the date of shipment, provided that the customer has stored and handled the product according to the information included on the product information sheet, website, and Certificate of Analysis. For living cultures, ATCC lists the media formulation and reagents that have been found to be effective for the product. While other unspecified media and reagents may also produce satisfactory results, a change in the ATCC and/or depositor-recommended protocols may affect the recovery, growth, and/or function of the product. If an alternative medium formulation or reagent is used, the ATCC warranty for viability is no longer valid. Except as expressly set forth herein, no other warranties of any kind are provided, express or implied, including, but not limited to, any implied warranties of merchantability, fitness for a particular purpose, manufacture according to cGMP standards,

typicality, safety, accuracy, and/or noninfringement.

Disclaimers

- This product is intended for laboratory research use only. It is not intended for any animal or human therapeutic use, any human or animal consumption, or any diagnostic use. Any proposed commercial use is prohibited without a license from ATCC.

While ATCC uses reasonable efforts to include accurate and up-to-date information on this product sheet, ATCC makes no warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. ATCC does not warrant that such information has been confirmed to be accurate or complete and the customer bears the sole responsibility of confirming the accuracy and completeness of any such information.

This product is sent on the condition that the customer is responsible for and assumes all risk and responsibility in connection with the receipt, handling, storage, disposal, and use of the ATCC product including without limitation taking all appropriate safety and handling precautions to minimize health or environmental risk. As a condition of receiving the material, the customer agrees that any activity undertaken with the ATCC product and any progeny or modifications will be conducted in compliance with all applicable laws, regulations, and guidelines. This product is provided 'AS IS' with no representations or warranties whatsoever except as expressly set forth herein and in no event shall ATCC, its parents, subsidiaries, directors, officers, agents, employees, assigns, successors, and affiliates be liable for indirect, special, incidental, or consequential damages of any kind in connection with or arising out of the customer's use of the product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, ATCC is not liable for damages arising from the misidentification or misrepresentation of such materials.

Please see the material transfer agreement (MTA) for further details regarding the use of this product. The MTA is available at www.atcc.org.

Copyright and Trademark Information

© ATCC 2023. All rights reserved.
ATCC is a registered trademark of the American Type Culture Collection.

Revision

This information on this document was last updated on 2023-04-15

iPSC-derived CD34+ Cells, BXS0117

ACS-7020

Product Sheet

Contact Information

ATCC
10801 University Boulevard
Manassas, VA 20110-2209
USA
US telephone: 800-638-6597
Worldwide telephone: +1-703-365-2700
Email: tech@atcc.org or contact your local distributor
