

ATCC Medium: 2174 Brain Heart Infusion Medium w/ 4mcg/ml Vancomycin

Base Medium

BHI Broth (BD 237500).....37.0 g
DI Water.....950 ml

Autoclave at 121°C.

Once base has cooled to 50°C, aseptically add the following:

Vancomycin Solution

Vancomycin.....0.004 g
DI Water.....50 ml

*If the reagent antibiotic powder indicates an assay value of potency, use that value to correctly calculate the amount of powder needed.

Filter sterilize.

*To make agar medium replace BHI Broth with 52 grams of BHI Agar (BD 211065)