

ATCC medium: 1861 ORI medium

Proteose Peptone No. 3 (BD 211693)	1.0 g
Yeast extract	1.0 g
Phytone Peptone (BD 211906)	0.5 g
Sodium thiosulfate	0.2 g
Sodium sulfite	0.05 g
Ferric citrate	0.4 g
Agar (for semi-solid medium)	3.0 g
Sea water	750.0 ml
Distilled water	250.0 ml

Adjust for final pH 7.6.
Autoclave at 121C for 15 minutes.