

ATCC medium: 2671 R2A (1:10 dilution) Normal Glucose Medium

Yeast extract.....	0.050 g
Proteose Peptone No. 3 (BD 211693).....	0.050 g
Casamino acids.....	0.050 g
Glucose.....	0.500 g
Soluble Starch.....	0.050 g
Sodium pyruvate.....	0.030 g
K ₂ HPO ₄	0.030 g
MgSO ₄ . 7H ₂ O.....	0.005 g
Agar.....	15.000 g
Distilled water.....	1.000 L

Adjust the final pH to 6.0. Autoclave at 121C for 15 minutes