

**ATCC medium: 80 Lactic bacteria broth**

Trypticase Peptone (BD 211921).....10.0 g  
Yeast autolysate.....5.0 g  
Sodium acetate.....12.0 g  
Glucose.....10.0 g  
Solution A (see below).....5.0 ml  
Solution B (see below).....5.0 ml  
Distilled water to.....1.0 L

Adjust pH to 5.1-5.3.  
Autoclave at 121C for 15 minutes.

*Solution A:*

$K_2HPO_4$  .....10.0 %  
 $KH_2PO_4$  .....10.0 %

*Solution B:*

$MgSO_4$  .....4.0 %  
 $NaCO_3$  .....0.2 %  
 $FeSO_4$  .....0.2 %  
 $MnSO_4$  .....0.2 %