

ATCC medium: 1834 Wilkins-Chalgren anaerobe medium (Oxoid CM643) with 2% yeast extract

Wilkins-Chalgren Anaerobe Broth (see below) supplemented with 2% (final concentration) yeast extract

Wilkins-Chalgren Anaerobe Broth:

Tryptone.....	10.0 g
Gelatin peptone.....	10.0 g
Yeast extract.....	5.0 g
Glucose.....	1.0 g
NaCl.....	5.0 g
L-Arginine.....	1.0 g
Sodium pyruvate.....	1.0 g
Menadione.....	0.5 mg
Hemin.....	5.0 mg
Distilled water.....	1.0 L

Adjust medium for final pH 7.1 (+/-) 0.2. Autoclave at 121C for 15 minutes.