

ATCC medium: 1394 *Halobacterium pharaonis* medium

Casamino Acids (BD 223050).....	15.0 g
Sodium citrate . 2H ₂ O.....	3.0 g
Glutamic acid.....	2.5 g
MgSO ₄ . 7H ₂ O.....	2.5 g
KCl.....	2.0 g
NaCl.....	250.0 g
Agar.....	20.0 g
Distilled water.....	1.0 L

Adjust pH prior to sterilization to 6.0. Autoclave at 121C for 15 minutes. Cool to 50C. Re-adjust pH aseptically to 8.5 after autoclaving.