

**ATCC Medium: 129 Nutrient Agar/Stabs/Broth with 0.5% NaCl**

Nutrient Agar (BD 213000).....23.0 g  
NaCl.....5.0 g  
DI Water.....1000 ml

Autoclave at 121°C.

**For Stabs**

Nutrient Broth (BD 234000).....8.0g  
NaCl.....5.0g  
Agar.....5.0g  
DI Water.....1000 ml

Autoclave at 121°C.

**For Broth**

Nutrient Broth (BD 234000).....8.0g  
NaCl.....5.0g  
DI Water.....1000 ml

Autoclave at 121°C.