

ATCC medium: 1999 *Azotobacter* medium (ATCC medium 240) with 1% sucrose

ATCC Medium 240 (see below) with 1% sucrose

ATCC Medium 240:

K ₂ HPO ₄	0.05 g
KH ₂ PO ₄	0.15 g
MgSO ₄ . 7H ₂ O	0.2 g
CaCl ₂	0.02 g
Na ₂ MoO ₄	0.002 g
FeCl ₃	0.002 g
Distilled water.....	1.0 L

For solid medium, add 1.5% agar.
Autoclave at 121C for 15 minutes.