

ATCC medium: 1536 *Anaerospirillum* medium

Polypeptone (BD 211910).....10.0 g
Yeast extract.....5.0 g
K₂HPO₄1.0 g
NaCl.....2.0 g
MgSO₄ . 7H₂O0.2 g
Na₂CO₃3.0 g
Glucose.....10.0 g
Distilled water.....1.0 L

Prepare and tube medium using anaerobic techniques and 100% CO₂ as gas phase. Autoclave at 121C for 15 minutes. Final pH of medium is 6.5.