

**ATCC medium: 1823 TSBY salt medium**

Trypticase Soy Broth (BD 211768).....30.0 g  
Yeast extract.....3.0 g  
NH<sub>4</sub>Cl.....0.25 g  
NaCl.....13.0 g  
MgCl<sub>2</sub> . 6H<sub>2</sub>O .....4.0 g  
MgSO<sub>4</sub> . 7H<sub>2</sub>O .....3.45 g  
KCl.....0.34 g  
CaCl<sub>2</sub> . 2H<sub>2</sub>O .....0.14 g  
Distilled water.....1.0 L

Adjust for final pH 7.2. Autoclave at 121C for 15 minutes.