

**ATCC medium: 291 Acidic tomato medium for *Leuconostoc***

Glucose.....10.0 g  
Peptone.....10.0 g  
Yeast extract.....5.0 g  
MgSO<sub>4</sub> . 7H<sub>2</sub>O .....0.20 g  
MnSO<sub>4</sub> . 4H<sub>2</sub>O .....0.05 g  
Tomato juice.....250.0 ml  
Agar (if needed).....15.0 g  
Distilled water.....750.0 ml

Adjust pH to 4.8.

Autoclave at 121C for 15 minutes.